Agar sangre/Agar Sal y Manitol/ Agar Mac Conkey

USO

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1. MacConkey, A. 1905. Lactose-fermenting bacteria in feces. *J. Hyg.* 5:333-379.
2. United States Pharmacopeial Convention, Inc. 1995. The United States pharmacopeia, 23 ed. The United States Pharmacopeial Convention, Rockville, M. D.
3. Food and Drug Administration. 1995 Bacteriological analytical manual, 8th ed. AOAC International. Gaithersburg, MD.
4. Gray, L.D. 1995 Escherichia, Salmonella, Shigella and Yersinia, p450-456. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
5. MacConkey J. H. 5:33. 1905. Joseph Md. State. Dept. Health. Procedures, 1960. European Pharmacopoeia. 4th Ed. 2002.
6. Song Y., Roumagna P., Weill F-X., Wain J., Dolecek C., Mazzone C.J., Holt K.E., y Achtman M. 2010 A multiplex single nucleotide polymorphism typing assay for detecting mutations that result in decreased fluoroquinolone susceptibility in Salmonella enteric serovars Typhi and Paratyphi A. *Journal of Antimicrobial Chemotherapy.* Vol. 65 1631-1641.
7. Jarvik T., Smillie C., Groisman E. A. y Ochman H. 2010 Short-Term Signatures of Evolutionary Change in the Salmonella enterica Serovar Typhimurium 14028 Genome. *Journal of Bacteriology.* Vol. 192 No. 2 560-567.
8. Porwollik S. 2011 Salmonella: From Genome to Function. Ed. Ilustrada Horizon Scientific Press. 300 Págs.
9. Farmacopea de los Estados Unidos Mexicanos: Suplemento para Dispositivos Médicos. 3a. Ed. -- México: Secretaría de Salud, Comisión Permanente de la Farmacopea de los Estados Unidos Mexicanos, 2014.
10. Leavitt, J.M. Naidorf and P. Shugaefsky. 1955. The undetected anaerobe in endodontic, a sensitive medium for detection of both aerobes and anaerobes. *The NYJ. Dentist.* 25:377-382
11. Curry, A. S., G.G., Joyce, and G.N. McEwen, Jr. 1993. CTFA Microbiology guidelines. The cosmetic, Toiletry and Fragrance Association. Inc. Washington, DC.
12. Brown, J.H. 1919. The use of blood agar for the study of streptococci. NY Monograph No. 9 The Rockefeller Institute for Medical Research.
13. The United States Pharmacopeia (USP XXIII) and The National Formulary (NF18). 1995 Sterility test, p. 1686-1690.
14. United States Pharmacopeia Convention Inc., Rockville, MD.
15. Vera and Power. 1980 In Lennette, Balows, Hausler and Tenover (ed.) Manual of clinical microbiology, 3rd ed. ASM, Washington, D.C.
16. Chapman, G.H. 1945. The significance of sodium chloride in studies of staphylococci. *J. Bacteriol.* 50:201-203.
17. Kloos, W.E., and T.L. Bannerman. 1995. Staphylococcus and Micrococcus. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C., Tenover, and R.H. Tenover (ed.) Manual of clinical microbiology, 6th ed. American Society for Microbiology. Washington, D.C.
18. United States Pharmacopeial Convention. 1995. The United States pharmacopeia, 2nd ed. The United States Pharmacopeial Convention, Rockville, M.D.
19. Hitchins, A.D., T.T. Tran, and J.E. MacCarron. 1995. Microbiology methods for cosmetics, p. 23.01-23.12. In Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, M.D.
20. McCulloch Am. j. Vet. Research, 8:173. 1947. Velilla, Faber, and Pelczar Am. J. Vet. Research, 8:275. 1947.
21. Jarvik T., Smillie C., Groisman E. A. y Ochman H. 2010 Short-Term Signatures of Evolutionary Change in the Salmonella enterica Serovar Typhimurium 14028 Genome. *Journal of Bacteriology.*
22. Rodríguez C.E. 2005 Bacteriología General: Principios y Prácticas de Laboratorio. Ed. Universidad de Costa Rica. 475 págs.

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